

Healthcare
furniture,
enhancing
lives.

Caremed Alrick.


Alrick EN7
Endless Bed Series

‘True floorline bed
with no towers and
endless possibilities’



The EN7 series bed is a variant of the EN9000 manufactured in Australia from imported components. It is the first ever universal bed that suits low care, high care and bariatric care, plus the innovative three-in-one bed. The EN7 uses Linak motors and has a low height of 100mm and has no towers!

Configuration Possibilities



Features

Electric hi/low	Auto regression backrest
Electric kneebend	75mm twin wheel castors
Electric forward & reverse Trendelenburg	Standard head and foot boards (fini blond laminate)

Specifications

Overall dimensions (mm)	2080 x 920 Single	2080 x 1090 King Single
Sleeping Surface dimensions (mm)	2000 x 900 Single	2000 x 1070 King Single
SWL	200kg *	
NET Bed Weight	110kg	
Low Position	100mm	
High Position	725mm	
Under bed Clearance	Clearview	
Angles (degrees)	Back	62
	Knee	35
	Reverse	15
	Trendelenburg	15
Warranty	Frame	10 yr
	Electrics	5 yr
	Handset & Castors	1 yr
Optional features	Underbed light	
	Battery backup	
	Dual braking	
	Upgrade to 300kg SWL *	
	End board options	

- 01 Electric forward & reverse Trendelenburg
- 02 Electric kneebend
- 03 Auto regression backrest
- 04 Standard head & footboard
- 05 10 year frame & 5 year electrical warranty
- 06 Electric Hi-low
- 07 75mm twin wheel castors

